



DIAGNOSTIC TESTING SUPPLIES

from MMS - A Medical Supply Company

Accustrip® Value+™ Strep A Rapid Test Strip

The Accustrip® Strep A Value+ Test Strip is a lateral flow, one-step immunoassay test for the rapid, qualitative detection of Group A Streptococcal antigen from throat swabs. The test is intended for use as an aid in the diagnosis of Group A Streptococcal infection. Test results in 5 minutes.

CLIA-Waived

JANID413M: 50 tests per box



Accutest® Integrated Strep A Test

The Accutest® Integrated Strep A rapid test device qualitatively detects the presence of Strep A antigen in throat swab specimens in 5 minutes. This rapid, easy to perform, 1-step test assists the physician in making treatment decisions during the patient's office visit.

CLIA-Waived

JANID440: 25 tests per box



Accutest® iFOBT Test

The Accutest® Single or Dual-Sample Immunological Fecal Occult Blood (iFOBT) Test is a rapid, qualitative, sandwich immunoassay for the determination of fecal occult blood in human stool. This innovative dual test allows two samples to be tested economically on the same test card. Testing more than one sample is recommended by the American Cancer Association.

CLIA-Waived

JANCS625: 25 Single-Well Test Devices and 25 collection tubes

JANCS601: 20 Single-Well Test Devices and 25 Mailers with collection material and tubes

JANCS612: 20 Dual-Well Test Devices and 25 Mailers with 2 collection materials and 2 tubes



800-473-2332 • www.mmsmedical.com

Your Clinical Laboratory Product Specialist



DIAGNOSTIC TESTING SUPPLIES

from MMS - A Medical Supply Company

Accustrip® Urinalysis Reader and Strips

With non-touch strip analysis and an easy-to-use interface, the Accustrip® URS Reader provides optimized patient data management. A simple touch screen coupled with great design and functionality make the Accustrip URS Reader a must-have for point-of-care urinalysis testing. Touch screen with graphic display, Instrument-read results, Automated, reliable and rapid results, Eliminate subjectivity of visual color interpretation, Reads 50 strips per hour, Patient ID or name entry with memory for 200 patient results, Automatic self-calibration. The Accustrip® urine reagent strips test for Blood, Urobilinogen, Bilirubin, Protein, Nitrite, Ketone, Glucose, pH, Specific Gravity and Leukocytes. Economical, Reliable accuracy. **CLIA-Waived** when read visually or with the Accustrip® URS Reader.

**JANUA870: 25 10-Parameter Urine Reagent Strips,
100 strips per bottle**

**JANUA850: URS Reader with Operator Manual and
110-120 volt AC power adaptor**

**JANUA812: URS Reader Starter Kit 12 – URS Reader and
12 bottles of 10-parameter Accustrip test strips**

**JANUA884: URS Reader Starter Kit 24 – URS Reader and
24 bottles of 10-parameter Accustrip test strips**

**JANUA832: Thermal Printer Paper for Accustrip URS Reader,
32 rolls per box**



Accutest® Urine Reagent Strips

Test for Glucose, Bilirubin, Ketone (acetoacetic acid), Specific Gravity, Blood, pH, Protein, Urobilinogen, Nitrite and Leukocytes in Urine. Accutest® URS-10 Urine Reagent Strips provide reliable sensitivity and are easy to read. **CLIA Waived**

JANUA710A: 100 strips per bottle



Accutest® Value+ One-Step hCG Test

The Accutest Value+ One-Step hCG Test is a rapid test for the detection of human chorionic gonadotropin (hCG) in urine specimens. This test is used to obtain a visual, qualitative result. Results in minutes. Easy-to-use cassette format. **CLIA-Waived**

JANPF535M: 40 tests per box



800-473-2332 • www.mmsmedical.com

Your Clinical Laboratory Product Specialist

MMS
A Medical Supply Company